




Global Health Crisis COVID-19

Surgical Face Masks:

Disposable flat type surgical face masks

- Type I
- Type II
- Type IIR

3 PROTECTION LAYERS



Performance requirements for medical facemasks (EN14683):

Test	Type I ^a	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm²)	< 29,4	< 29,4	< 49,0
Splash resistance pressure (kPa)	Not required	Not required	≥ 16,0
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

a. Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

Features:

- **CE marked** face masks
- Easy-wear ear hanging type
- Great fit
- Three protection layers

Manufacturer:

CHIMIN HEALTH MANAGEMENT, Co. Ltd., Industrial group of the healthcare sector to which LINEAR belongs, manufactures surgical masks in one of its factories in the city of Taizhou, with a **daily production capacity of 400,000 units.**

