



**Thermo**  
SCIENTIFIC  
Authorized Distributor

## Heracell VIOS 160i CO2 Incubator



### Features:

- 165 litre capacity
- Decontamination cycle 180°C on all surfaces
- THRIVE active airflow system with HEPA filter-ISO CLASS 5
- 2-year warranty
- Air jacketed for faster recovery times
- Advanced in-situ sensor technology
- Electro-polished surface to minimize bacterial growth
- Interactive touchscreen control panel i-CAN
- Temperature range ambient +3°C to 55°C
- CO2 control from 1-20%
- Integral humidity system with water level sensor
- Data capturing with USB output




**Part # 51030287**

A CO2 incubator specifically designed to deliver the performance, reliability, ease of operation and value required to support a range of culturing needs from basic research to demanding leading-edge applications.



**Offer Valid Till End of August 2020 - Contact Us Now!**



Separation Scientific SA (Pty) Ltd

Tel : 011 794 4430  
Email : [info@sepsci.co.za](mailto:info@sepsci.co.za)  
Web: [www.sepsci.co.za](http://www.sepsci.co.za)