



Thermo
SCIENTIFIC
Authorized Distributor



Part # 50129845PROMO

Thermo Smart2Pure12 UV/UF with Tank


- Dual Type I and Type II system
- 18.2MΩ purity 0.055uS/cm
- Bacteria count < 1cfu/ml
- Particle content <1
- TOC 1-5ppb (with UV)
- Endotoxins 0.001EU/ml (with Ultra filter)
- Permeate production at 15°C – 12 liter per hour



**Delivery and Installation in Gauteng included.
For all other regions - POR**



- Storage tank 60L with conical bottom
- Recirculation between unit/tank to maintain purity
- Microprocessor control with filter limit alarms
- Water qualities to meet/exceed ASTM I and II, CAP, ISO3696, BS3997 and CLSI standard
- 185/254 nm UV lamp
- RS232 interface for GLP printouts
- Incorporates 5 stage purification platform with Prefiltration, RO, Deionization, UV and UF in one complete unit
- Cost effective maintenance and service by user



Offer Valid Till End of September 2020 or while stocks last - Contact Us Now!

Separation Scientific SA (Pty) Ltd

Tel : 011 794 4430
Email : info@sepscico.co.za
Web: www.sepscico.co.za