



SCICOOOL 1140R Glass Door Medical Cooler with Low E Glass

Thermo
SCIENTIFIC
Authorized Distributor



Part # SCICOOOL1140R

Features

- LED Lights: Standard
- Electronic Controller (Thermostat with Display): Standard
- Over/Under Voltage Compressor Protection: Standard
- Lock: Standard
- 12 Month Warranty

Dimensions


- External: (H) 2020mm x (W) 1140mm x (D) 635mm
- Internal (Nett): (H) 1380mm x (W) 1050mm x (D) 550mm



Delivery and Installation in Gauteng included.
For all other regions – Price On Request

Specifications

- Gross Capacity (L): 890
- Nett Capacity (L): 750
- Shelves: 8 (extra shelves optional)
- Refrigerant: R134a
- Temp Range: 2°C to 8°C (Factory preset to 4°C)
- Max Power Consumption (watt): 637
- Voltage/Hertz: 230/50
- Compressor Size (hp/cc): 3/8
- Weight Empty (kg): 180
- Swinging Door



Separation Scientific SA (Pty) Ltd

Tel : 011 794 4430
Email : info@sepsci.co.za
Web: www.sepsci.co.za

Offer Valid While Stocks Last - Contact Us Now!