



MODEL ULT25NEU

PORTABLE ULTRA-LOW TEMPERATURE STORAGE

Safeguards temperature-sensitive medical treatments regardless of patient location.

- ✔ Lightweight, portable design makes this a perfect appliance for clinical point-of-use cold storage
- ✔ Superior Stirling engine reliability scales down to a remarkably small footprint and weight (21 kg, 46 lbs), without sacrificing performance
- ✔ Provides industry's widest ultra-low temperature range of -20°C to -86°C
- ✔ Adaptable power input provides operational flexibility between standard AC main outlets, worldwide, as well as vehicle-based 12V DC power or external battery devices
- ✔ Utilizes 100% natural refrigerants to support sustainable, environmentally friendly operations
- ✔ Optional Sense Anywhere wireless temperature logging allows for remote monitoring via cloud-based software
- ✔ Ample storage capacity for thousands of 2mL vials within a surprising compact frame



STIRLING
POWER

Industry's widest ULT temperature range -20°C to -86°C

Uniquely portable and reliable ultra-low ULT25NEU storage benefits scientists, clinicians, and patients alike.

Contact us now for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za