




Scientific Cold Storage Solutions



SCI-COOL 1140F

Glass Door Medical Freezer with Low E Glass



<i>Features</i>	SCI-COOL 1140F
Automatic Frost-free System Design	Standard
Heated Glass Door	Standard
LED Lights	Standard
Electronic Controller (Thermostat With Display)	Standard
Over/Under Voltage Compressor Protection	Standard
Lock	Standard
Front Ventilation	Standard
Hinge Doors	2
Available in Glass or Solid Doors	Yes
<i>Specifications</i>	
Gross Capacity	879 Ltr
Nett Capacity	739 Ltr
Shelves	10 (extra shelves optional)
Refrigerant	R 404
<i>Temp Range</i>	-15°C to -25°C
Voltage / Hertz	220/50
Weight Empty (kg)	195
<i>Dimensions</i>	
External	(H) 2020mm x (W) 1140mm x (D) 635mm
Internal (Nett)	(H) 1380mm x (W) 1050mm x (D) 550mm

**TINTED
DOORS
AVAILABLE**

Contact us for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za