




Scientific Cold Storage Solutions



SCI-COOL 580F

Glass Door Medical Freezer with Low E Glass



<i>Features</i>	SCI-COOL 580F
Automatic Frost-free System Design	Standard
Heated Glass Door	Standard
LED Lights	Standard
Electronic Controller (Thermostat With Display)	Standard
Over/Under Voltage Compressor Protection	Standard
Lock	Standard
Front Ventilation	Standard
<i>Specifications</i>	
Gross Capacity	251 Ltr
Nett Capacity	207 Ltr
Shelves	3 (extra shelves optional)
Refrigerant	R 290
<i>Temp Range</i>	
	-15°C to -25°C
Voltage / Hertz	220/50
Weight Empty (kg)	90
<i>Dimensions</i>	
External	(H) 1350mm x (W) 580mm x (D) 635mm
Internal (Nett)	(H) 840mm x (W) 490mm x (D) 510mm

**TINTED
DOORS
AVAILABLE**

Contact us for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za