



avantor™
delivered by **VWR™**



Benchtop pH/MV meter

Features automatic or manual temperature compensation, automatic or manual ranging, automatic buffer recognition, digital and analog outputs, min/max and 99 data points, data hold, 5 calibration points, and electrode status indication. Ordering information: 5-year warranty.

Check this out!

Sporicidin® disinfectant wipes and towelettes are EPA registered disinfectant wipes recommended for use in hospitals, and medical, dental, and veterinary clinics. Sporicidin® Disinfectant towelettes are used for hospital disinfecting, compliant for medical instrument disinfecting prior to sterilization

**Sporicidin®
Disinfectant Wipes**



Spectrophotometers

- 1. VWR® Double Beam UV-Vis Spectrophotometer**
- 2. VWR® UV-Vis Scanning Spectrophotometer**

Supplied with 4 optical glass cells, 2 quartz cells (excludes model V3000-PC), dust cover, USB cable, USB (8 GB) stick, UV-Vis analyst software for PC control, and 4-position cell holder (10 mm).

PerfeCTa® qPCR ToughMix®, Low ROX™, QuantaBio

PerfeCTa qPCR ToughMix, Low ROX is a 2X concentrated ready-to-use reaction cocktail for PCR amplification of DNA templates that overcomes many known inhibitors of PCR often present in crude samples extracted from environmental specimens, plant tissues, or animal tissues.



Separation Scientific

Tel: 011 794 4117 | Email: Info@sepscico.za | Web: www.sepscico.za



avantor™
delivered by **VWR™**



VWR® Advanced Dry Block Heaters (1, 2 or 4 blocks)

500 nm, 400 mg/dL, Trinder
Methodology –
Colorimetric
Reagents are designed for
maximum stability, ease-of-
use and are optimized for
today's chemistry analyzers.
Substrates are ready to use
liquid reagents.

Multi-purpose units are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points, and a variety of other laboratory procedures

RoHS certified. TUV listed. Include a Troemner ISO/IEC accredited calibration certificate traceable to NIST. The ISO/IEC 17025 accredited calibration certificate includes multiple data points within the temperature and time function with the associated uncertainties. Troemner is accredited by NVLAP under laboratory code 105013-0 to perform ISO/IEC 17025 calibrations.



Glucose LiquiColor® Test (Enzymatic), 1L



Contact **Sep Sci** for more information on our Best Seller range.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za