



Quantitative Standards Controls



NATtrol Adenovirus Type 1 Stock is formulated with purified, intact viral particles that have been chemically modified to render them non-infectious and refrigerator stable. Each vial contains 1.0 mL of NATtrol Adenovirus Type 1 stock in a proprietary matrix at the target concentration.

NEW

QUANTITATIVE MS2 STOCK

*NATtrol MS2 Bacteriophage
stock is now available!*
NATMS2-STQ - 1.0 mL

NATtrol Adenovirus Type 1 stock is designed to evaluate the performance of nucleic acid tests for determination of the presence of Adenovirus Type 1 nucleic acid. NATtrol Adenovirus Type 1 stock enables laboratories to monitor test variation, lot-to-lot test kit performance, operator variation, and can provide assistance in identifying random or systemic error.

Real-Time PCR



Adenovirus

NATTROL

Positive

Contact us now for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za