




## NEW FLU CONTROL

NATtrol Flu/RSV/SARS-CoV-2 External Run Control



NATtrol Flu/RSV/SARS-CoV-2 External Run Control (NATFRC-ERC) is formulated with purified, intact viral particles that have been chemically modified to render them non-infectious and refrigerator stable. Each control pack contains 6 x 0.5 mL of NATtrol inactivated virus.

NATtrol Flu/RSV/SARS-CoV-2 External Run Control is designed to evaluate the performance of nucleic acid tests for determination of the presence of respiratory virus nucleic acids. NATFRC-ERC enables laboratories to monitor test variation, lot-to-lot test kit performance, operator variation, and can provide assistance in identifying random or systemic error.

### NATTROL

#### NATFRC-ERC

contains 6 x 0.5 mL of NATtrol inactivated virus. Each vial has a combination of Flu A / Flu B / RSV B / Sars-CoV-2. Target concentration levels are 15,000 cp / mL

Platform : GeneXpert  
Test: SARS-COV-2 / +2° to +8° C.

Contact **Sep Sci**  
for more information.

**Separation Scientific**

Tel: 011 794 4117 | Email: [Info@sepsci.co.za](mailto:Info@sepsci.co.za) | Web: [www.sepsci.co.za](http://www.sepsci.co.za)