




MONKEY POX EXTERNAL RUN CONTROL



MONKEYPOX EXTERNAL RUN CONTROL IS A FULL PROCESS, WHOLE ORGANISM, INACTIVATED CONTROL

- ✓ Refrigerator stable
- ✓ Ready to use
- ✓ Mimics a true clinical specimen

NATMPXV-ERC CONTAINS 6 X 0.5 ML VIALS

NATtrol Monkeypox Virus External Run Control is formulated with purified, intact viral particles that have been chemically modified to render them non-infectious and refrigerator stable. Each control pack contains 6 x 0.5 mL of NATtrol Monkeypox virus external run control in a proprietary matrix at the target concentration listed.

NATtrol Monkeypox virus external run control is designed to evaluate the performance of nucleic acid tests for determination of the presence of Monkeypox virus DNA. NATtrol Monkeypox virus external run control enables laboratories to monitor test variation, lot-to-lot test kit performance, operator variation, and can provide assistance in identifying random or systemic error.

NEW

CONTACT US NOW

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za