



Gastrointestinal Diagnostic Solutions

Gastrointestinal diagnostics have long been a pillar of our company's mission. From H. pylori to C. difficile, to Foodborne, our solutions provide diagnostic accuracy you expect with simplified efficiencies to optimize workflow. We provide flexible testing solutions to fit your unique situation and variable testing demands.

Premier®

The Premier testing platform is a microwell enzyme immunoassay (ELISA) technology that delivers optimal performance and accurate results. Premier® products use detachable microwells.

Platform Specifications:

1. Microwell plates allows up to 96 determinations
2. ELISA technology lends itself for high volume laboratories
3. Can be read visually and on automated systems
4. Fast and simple detection with streamlined test procedures
5. Refrigerated kit storage

• Premier Adenoclone	696007
• Premier® Adenoclone® 40/41	696006
• Premier® CAMPY	618096
• Premier® C.difficile GDH	611096EU
• Premier® Coccidioides	603096
• Premier® EHEC	602096
• Premier® H. pylori	608096
• Premier® Platinum HpSA® PLUS	606096
• Premier® Rotaclone®	601396
• Premier® Toxins A & B	696004
	616096

Healthcare Acquired Infections

- C. difficile

Gastrointestinal

- Calprotectin
- H. pylori
- Rotavirus
- E.coli
- Adenovirus
- Campylobacter

Parasitology

- Cryptosporidium/Giardia
- Para Pak

Curian

- Curian HpSA Assay 760130

Immunocard®

- Immunocard C.difficile GDH 716050
- Immunocard® H. pylori 710030
- Immunocard® Mycoplasma 709030
- Immunocard® Toxins A&B 712050

Immunocard STAT!®

Immunocard STAT! is a rapid, lateral flow immunoassay technology that delivers fast results with optimal performance and confidence.

• Immunocard STAT! C. difficile GDH-AB	750520
• Immunocard STAT!® C. difficile GDH-AB External Controls	750501
• Immunocard STAT!® CAMPY	751530
• Immunocard STAT!® CGE	751420
• Immunocard STAT!® CGE Positive Control	751403
• Immunocard STAT!® Crypto/Giardia	750830
• Immunocard STAT!® Crypto/Giardia Positive Control	750820
• Immunocard STAT!® E. coli O157 PLUS	750530
• Immunocard STAT!® EHEC	751630
• Immunocard STAT!® HpSA® (CLIA Waived)	750220
• Immunocard STAT!® Rotavirus	750030
• Immunocard STAT!® Rota-Adeno-Noro2	750320
• Immunocard STAT!® Rota-Adeno-Noro2 External Controls	750301



We want to help providers make better **diagnostic decisions** with a focus on **gastrointestinal conditions.**

Contact us now for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za