



NATTROL

● RESPIRATORY VERIFICATION PANEL

Controls and Standards

**NATtrol Respiratory
Verification Panel** (qualitative)
is formulated with purified,
intact bacterial cells and viral
particles.

The microorganisms have
been chemically modified to
render them non-infectious
and refrigerator stable.
NATRSP-BIO contains
23 x 0.6 mL vials of bacterial
and viral NATtrol and
5 x 1.8 mL vials of negative
(matrix only).



- **NATRSP-BIO**
(23 x 0.6 mL, 5 x 1.8 mL)
- **NATRSPC-BIO**
(12 x 0.3mL)

External Run Controls

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za