

NEW NATrol respiratory verification panel mini



NATtrol Respiratory Verification Panel Mini* (qualitative) is formulated with purified, intact viral particles that have been chemically modified to render them non-infectious and refrigerator stable. NATRMP-BIO contains 6 x 0.6 mL vials of viral NATtrol and 2 x 1.8 mL vials of negative (matrix only) as listed in Table 1. The panel members are supplied in a proprietary matrix.

NATtrol Respiratory Verification Panel Mini is designed to evaluate the performance of nucleic acid tests for determination of the presence of viral nucleic acids (from organisms listed in Table 1). NATtrol Respiratory Verification Panel Mini can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.

Platform: BioFire Spotfire; Test: Respiratory / +2° to +8° C.

- ✓ Quality Control Panels
- ✓ Validation-Verification
- ✓ NATtrol Respiratory Verification Panel Mini

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za