



Quality Control Panels Validation-Verification

NATtrol Global Fever Panel

NATtrol Global Fever panel NATGFP-BIO is now available

This panel (22 x 0.3 mL) has been tested with the BioFire® Global Fever Panel assay and provides all expected results.

NATtrol Global Fever Panel is formulated with purified, intact bacterial cells and viral particles. The microorganisms have been chemically modified to render them non-infectious and refrigerator stable. **NATGFP-BIO** contains 22 x 0.3 mL vials of bacterial and vira NATtrol. The panel members are supplied in a proprietary matrix.



NATtrol Global Fever Panel is designed to evaluate the performance of nucleic acid tests for determination of the presence of bacterial, protozoan and viral nucleic acids. **NATtrol Global Fever Panel** can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel (+2° to +8° C).

Chikungunya | Dengue | E. coli | NATTROL | Positive | West Nile Virus | Yellow Fever

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za