



GASTRO




INTESTINAL CONTROLS

The nervous system that regulates our gut is often called **the body's second brain.**

NATGIP2-QIA
NATtrol GI Verification panel

Quality Control Panels



Make the smart choice when selecting GI controls and panels. **NATtrols** can be used across various molecular testing platforms and assays

NATtrol GI (Gastrointestinal) Verification Panel 2 is formulated with purified, intact cells (bacterial, fungal or parasites) and viral particles. The cells and viral particles have been chemically modified to render them non-infectious and refrigerator stable. NATGIP2-QIA contains 22 x 0.25 mL vials of bacterial, fungal, parasite and viral NATtrol and 4 x 1.0 mL of negative control (matrix only). The panel members are supplied in a proprietary matrix.

NATtrol GI Verification Panel 2 is designed to evaluate the performance of nucleic acid tests for determination of the presence of bacterial, parasitic and viral nucleic acids. NATtrol GI Verification Panel 2 can also be used for validation of clinical assays, development of diagnostic tests, and training of laboratory personnel.

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za