



NATtrol RSP Sore Throat Supplement

NATtrol RSP (Respiratory Verification Panel Sore Throat Supplement) is now available!

NATRSPX-BIO (Qualitative)
contains 2 x 0.6 mL
vials of bacterial
NATtrol

NATRSPX-BIO

NATtrol RSP (Respiratory Verification Panel) Sore Throat Supplement (qualitative) is formulated with purified, intact bacterial cells. The microorganisms have been chemically modified to render them non-infectious and refrigerator stable. NATRSPXBIO contains 2 x 0.6 mL vials of bacterial NATtrol. The panel members are supplied in a proprietary matrix.

NATtrol RSP Sore Throat Supplement is designed to evaluate the performance of nucleic acid tests for determination of the presence of bacterial nucleic acids from *S. dysgalactiae* and *S. pyogenes*. NATtrol RSP Sore Throat Supplement can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.



Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za