

Finnpipette™ F1 Variable Volume Pipettes

Enjoy the enhanced volume adjustment on Thermo Scientific™ Finnpipette™ F1 Variable Volume Pipettes, which is now silent and effortless with textured adjustment knob for a firmer feel when setting volumes. The adjustable finger rest and lightweight design maximize comfort without sacrificing performance.

Product Code	Description	Qty
4641010N	* Finnpipette F1 0.2-2µl, micro	6
4641020N	* Finnpipette F1 0.5-5µl, micro	4
4641030N	* Finnpipette F1 1-10µl, micro	6
4641050N	* Finnpipette F1 2-20µl, micro	4
4641070N	Finnpipette F1 10-100µl	4
4641080N	Finnpipette F1 20-200µl	1



* Micro range pipettes are only compatible for use with the Finntip micro brand of tips

Finnpipette™ F2 Variable Volume Pipettes

Meet the challenges of highly demanding applications where robustness and durability are essential with the sleek, rugged design of Thermo Scientific™ Finnpipette™ F2 Variable Volume Pipettes. Built with tough PVDF components that stand up to harsh chemicals and the damaging effects of UV light, these variable volume pipettes are fully autoclavable without disassembling. They feature a wide selection of expressly designed tips for optimal performance, precision and accuracy.

Product Code	Description	Qty
4642010	* Finnpipette Fp F2 0,2-2 µl µTip	5
4642030	* Finnpipette Fp F2 1-10 µl µTip	7
4642050	* Finnpipette Fp F2 2-20 µl µTip	3
4642070	Finnpipette Fp F2 10-100 µl	4
4642080	Finnpipette Fp F2 20-200 µl	5
4642090	Finnpipette Fp F2 100-1000 µl	3



* Micro range pipettes are only compatible for use with the Finntip micro brand of tips

Offer valid while stock lasts!!

ORDER HERE

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za

Offer valid while stock lasts!!

Finnpipette™ Novus Electronic Pipettes

Thermo Scientific Finnpipette Novus single and multichannel electronic pipettes. Finnpipette® Novus delivers performance - with greater ease-of-use and less risk of injury. Unique to Novus, the adjustable index-finger trigger action prevents common repetitive stress injuries caused by thumb-driven pipettes, and the easy-to-follow prompts guide Novus users through 10 pipetting functions. For a full range of applications, Novus will be your new favorite pipette.

Product Code	Description	Qty
46200400	Finnpipette FP Novus 10-100µl	1
46200500	Finnpipette FP Novus 30-300µl	1
46200600	Finnpipette FP Novus 100-1000µl	1
46300300	Finnpipette Novus MCP 12 5-50µl	4
46300500	Finnpipette Novus MCP 12 30-300µl	1

Product Code	Description	Qty
9400230	Finntip 250 Universal, 1000/bag	30
9400303	Finntip 10, 10x96/rack, sterile	3
9400310	Finntip 10, 1000/bag	10
94052000	Finntip Filter 10ul Micro, 10x96/rack, s	10
94052150	Finntip Filter 20ul, 10x96/rack, sterile	10
94052200	Finntip Filter 100ul, 10x96/rack, sterile	10
94052300	Finntip Filter 200ul, 10x96/rack, sterile	14
94052600	Finntip Filter 10ml, 5x24/rack, sterile	5
94060810	Finntip Flex 1200, 10x96/rack	3
94060820	Finntip Flex 1200, 1000/bag	3
94300120	Finnpipette Labtip Yellow, 1000/bag	17
94300220	Finnpipette Labtip Blue, 1000/bag	39
9420340	Finnpipette rotating stand	4
95128093	Reagent Reservoir 25 ml standard 100 pcs	9



Pipette Tips



ORDER HERE

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za

F1-ClipTip™ Variable Volume Single Channel Pipettes

Ensure a secure fit with ClipTip technology, which locks tips firmly to the pipette with a complete seal so they will not loosen or fall off. Lightweight handle design combined with minimal tip attachment and ejection forces will increase the comfort in pipetting. Feel the difference.

Product Code	Description	Qty
4641180N	F1-ClipTip 2-20µl	3
4641200N	F1-ClipTip 10-100µl	3
4641210N	F1-ClipTip 20-200µl	2
4641220N	F1-ClipTip 30-300µl	1
4641230N	F1-ClipTip 100-1000µl	2
4641310N	F1-ClipTip 0.1-2 µl	3
94420213	ClipTip 20 Filter, 10x96/rack, Sterile	5
94420313	ClipTip 200 Filter, 10x96/rack, Sterile	1
94420713	ClipTip 1000 Filter, 8x96/rack, Sterile	4



Offer valid while stock lasts!!

ORDER HERE

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za