



# BIOFIRE GLOBAL FEVER PANEL



NATtrol™ Global Fever Panel is formulated with purified, intact bacterial cells and viral particles

NATtrol Global Fever Panel is formulated with purified, intact bacterial cells and viral particles. The microorganisms have been chemically modified to render them non-infectious and refrigerator stable. NATGFP-BIO contains 22 x 0.3 mL vials of bacterial and viral NATtrol. The panel members are supplied in a proprietary matrix.

NATtrol Global Fever Panel is designed to evaluate the performance of nucleic acid tests for determination of the presence of bacterial, protozoan and viral nucleic acids. NATtrol Global Fever Panel can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.

## NATGFP-BIO

NATGFP-BIO contains 22 x 0.3 mL vials of bacterial and viral NATtrol

ORGANISMS: Chikungunya, Dengue, E. coli, West Nile Virus, Yellow Fever



Quality  
Control  
Panels

**Separation Scientific**

Tel: 011 794 4117 | Email: [Info@sepsci.co.za](mailto:Info@sepsci.co.za) | Web: [www.sepsci.co.za](http://www.sepsci.co.za)