

x
x
x



thermo
scientific

Authorized Distributor

GET READY FOR THIS

**STOCK
CLEARANCE
SPECIAL**

while stocks last

x
x
x



Pico 21

**For micro-volume protocols
such as nucleic acid, protein lysate
preparation and PCR reaction set-up.
All in a remarkably small footprint!**

Thermo Scientific™ Pico 21 instrument complete unit with:

- 24 x 1.5/2ml micro rotor with ClickSeal biocontainment lid
- Unique Dual row 18 x 2ml/0.5ml rotor with screw on lid
- Max. RCF 21 000 x g
- Max. Speed 14,800rpm

**R 54,190.00
excl vat**

- **Net Weight: 10.5kg**
- **Capacity: 24 x 1.5/2mL**
- **Noise Level: 56dbA**
- **Warranty: 2 years**

- **Electrical Requirements: 230V 50/60Hz**
- **Run Time: 1 to 99min (1 min increments)**
- **Standards: IEC 61010-1, IEC 61010-2-020**
- **Dimensions (HWD): 22.5 x 24.3 x 35.2cm**

Features

ORDER HERE

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za

x
x
x