


Quality Control Panels



NATRSP-DIA

NATtrol Respiratory Verification Panel can be used to develop procedures for verification of performance of molecular assays that detect the presence of bacterial and viral nucleic acids from organisms. The panel can also be used for training and evaluation of operator proficiency.

NATRSPC-DIA

NATtrol Respiratory Panel Control is a qualitative external control intended to be used with molecular assays that detect nucleic acids from organisms. The control can be used to monitor the steps of extraction, amplification, and detection, and also for training and evaluation of operator variation.

NATTROL

**Two NEW
qualitative respiratory
panels for you,
calibrated for the
Diasorin liaison
plex platform**

1

**NATRSP-DIA Respiratory
Verification Panel**

2

**NATRSPC-DIA Respiratory
Panel Control**



Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za