



Super Sale

Medifuge



LIMITED STOCK

Ideal for: Clinical and diagnostic laboratories, wards, veterinary labs and small hospitals. Cell culture pelleting, regenerative therapy.

- Innovative Dual-Spin™ rotor for fixed angle & swingout applications.
- **8 x 15ml (fixed angle) and 8 x 10ml (swing out).**
- 300–4900rpm (100rpm increments).
- Time setting 1–99minutes, HOLD.
- Programmable with 4 quick access buttons.
- Bright LED display with microprocessor control.
- Acceleration and deceleration options.
- Maintenance free induction motor.
- **UL, CE, RoHS, CSA, IVD, FDA class 1.**

2 YEARS WARRANTY

R 29,995.00

Cat No: **75008800** Offer includes

- Dual-Spin rotor
- Fixed angle buckets (set8)
- Swinging buckets (set8)
- Green short spacers (set8)
- Yellow long spacers (set8)

LIMITED STOCK

Cat No: **75004081** Offer includes:

R 17,780.00

- 12 x 1.5/2ml rotor
- 32 x 0.2ml tube rotor (single or 4x8 strips)
- 0.2ml tube adaptor (set of 12)
- 0.5ml tube adaptor (set of 12)



2 YEARS WARRANTY

MySpin 6



R 7,255.00

Cat No: **75004061** Offer includes

- 6 x 1.5/2ml rotor
- 16 x 0.2ml tube rotor (singled or 2x8 strips)
- 0.2ml tube adaptor (set of 6)
- 0.5ml tube adaptor (set of 6)

LIMITED STOCK

2 YEARS WARRANTY

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za