





Ready to use media

DMEM 20% FBS

DMEM w/ 4.5 g/l Glucose, w/ L-Glutamine, w/ Sodium Bicarbonate, 3.7 g/l NaHCO₃ 500ml

0.2% PenStrep Fungizone
0.2% Normase

Available Now

Prepared By:	Date prepared: 20-02-2024
Lot Number: C465887	Expiry date: 20-02-2025
Storage: 2-8 C	

**Save Time
in the Lab
and Save
Money!!**

**READY
TO USE
MEDIA**

Cellonex now supplies media that is made up for you and ready to use!!

All common **Cell Culture Media** are offered as ready to use, e.g. DMEM, RPMI 1640, and DMEM:HAM's F12. Prepared with the concentration of FBS that you need already added. (2%, 5%, 10%, 20%). *Premium grade serum is used.



♥ PenStrep/Fungizone and Anti-mycoplasma reagents are added(0.2%). This is the media and serum that **Cellonex** uses to grow up their cells in the repository. Our cells love it. ♥



**Premium Grade
Foetal Bovine
Serum**

We also offer **FBS aliquots** (Premium grade) in 100 ml 250ml and 500ml - **No permits required.** Only buy what you need and save your grant money for other essential reagents. All media is prepared in sterile conditions in a clean room where strict GLP is implemented.