

VDRL ANTIGEN MR

Determination of plasma reagins


SLIDE TEST

VDRL MR 250 tests:
Kit does not include controls,
however Positive and Negative
controls can be purchased
separately if required

**SYPHILIS**

The assay is a modified VDRL Test. The **VDRL Antigen MR** is a non-treponemal preparation specially developed for the rapid detection and semi-quantitation by flocculation on a slide of plasma reagins, a group of antibodies directed against tissue components produced by almost every patient infected with *T. pallidum*. The assay is performed by testing the antigen –an association of lipid complexes and choline chloride– against not inactivated samples.

- ✓ Syphilis is caused by infection with the bacterium *Treponema pallidum* which can be transmitted congenitally or by sexual contact.
- ✓ The VDRL MR test permits a rapid screening of large numbers of persons so that reactors can be given treatment.
- ✓ VDRL MR test has a high diagnostic value on a tentative diagnosis made on the basis of case history and clinical findings.



For in vitro diagnostic use only

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za

