



ZeptoMetrix®

Accuracy Matters

H5N1 IS HERE. YOUR CONTROLS SHOULD BE TOO!



Product # 0831198

NATtrol Influenza A H5N1 Control is:

Non-infectious

Ready-to-use

Target 1,000,000 gene copies/mL

Designed to support accurate assay performance



Influenza A H5N1 continues to present zoonotic and pandemic threat. This virus has demonstrated the ability to cause severe disease in humans with high case fatality rates, and has shown expanded host range and geographic spread.

ZeptoMetrix offers laboratory a comprehensive **H5N1 quality control material** that requires RNA extraction, unlike traditional naked RNA products. This novel solution includes **three encapsulated H5N1 gene segments** within a phage-like particle, ensuring realistic validation of extraction, amplification, and detection steps.

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za