




Antigens  
Vials -ve

# Treponema Pallidum Antigen

Looking For  
A Good  
Quality  
Treponema

Treponema pallidum Antigen (Nichols Strain). Cells of Nichols Strain Treponema pallidum organisms supplied in phosphate buffered saline with 0.1% sodium azide and 15 % glycerol. Appearance: white opalescent suspension, frozen liquid. Source: Rabbit testis. Suitable for use in **Serodiagnostic Testing** used for the diagnosis of Syphilis.



Contact **Sep Sci** for more info.

**Separation Scientific**

Tel: 011 794 4117 | Email: [Info@sepsci.co.za](mailto:Info@sepsci.co.za) | Web: [www.sepsci.co.za](http://www.sepsci.co.za)