



LINEAR



GASTROENTEROLOGY



KROMA iT

- i-FOB Turbidimetric Fecal Occult Blood
- i-TRA Turbidimetric Transferrin
- i-CPF Turbidimetric Calprotectin



- 150 tests/hour, up to 250 tests/hour with ISE module (optional)
- 9 reagent positions
- 27 sample positions
- Compatible with Mono, Double, and Triple reagent methodologies
- Cap Piercing

Automatic Immunological Analyzer



KROMA Plus iT

- i-FOB Turbidimetric Fecal Occult Blood
- i-TRA Turbidimetric Transferrin
- i-CPF Turbidimetric Calprotectin



- 240 tests/hour, up to 360 tests/hour with ISE module (optional)
- 35 reagent positions
- 49 sample positions
- Compatible with Mono, Double, and Triple reagent methodologies
- Cap Piercing.

Automatic Immunological Analyzer



LIDA 300

- i-FOB Turbidimetric Fecal Occult Blood
- i-TRA Turbidimetric Transferrin
- i-CPF Turbidimetric Calprotectin



- 420 tests per hour
- Clinical chemistry analyzer for lab tests on bodily fluids (glucose, cholesterol, enzymes, etc.)

Clinical Chemistry Analyzer

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za