

x
x
x



thermo
scientific

Authorized Distributor

GET READY FOR THIS

STOCK CLEARANCE SPECIAL

x
x
x

Microcentrifuge Sale

x
x
x

For micro-volume protocols
such as nucleic acid, protein lysate
preparation and PCR reaction set-up.
All in a remarkably small footprint!



**75002410
PROMO**

Thermo Scientific™ Pico 17 instrument:

- Complete with 24 x 1.5/2ml rotor
- Biocontainment certified ClickSeal lid
- RCF 17 000 rcf xg at 13 300 rpm
- 2 year warranty with Thermo Scientific leading performance and support reputation.

- Capacity: 24 x 1.5/2mL
- Noise Level: 56dbA
- Warranty: 2 years

- Electrical Requirements: 230V 50/60Hz
- Run Time: 1 to 99min (1 min increments)
- Standards: IEC 61010-1, IEC 61010-2-020

Features

ORDER HERE

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za

x
x
x